# For Research Use Only. Not for use in diagnostic procedures.



#### MONOCLONAL ANTIBODY

# Alexa Fluor® 488 labeled Mouse IgG2a isotype control

Code No. Clone Subclass Quantity Concentration M076-A48 6H3 Mouse IgG2a  $\kappa$  100  $\mu L$  1 mg/mL

**SOURCE:** This antibody was purified from hybridoma (clone 6H3) supernatant using protein A agarose. This hybridoma was established by fusion of mouse myeloma cell P3U1 with Balb/c mouse lymph nodes immunized with KLH.

**FORMULATION:** 100 μg IgG in 100 μL volume of PBS containing 1% BSA and 0.09% NaN<sub>3</sub>.

\*Azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush plenty of water when disposing materials containing azide into drain.

**STORAGE:** This antibody solution is stable for one year from the date of purchase when stored at 4°C.

**REACTIVITY:** No specific binding is detected on human peripheral blood lymphocyte, monocyte and granulocyte.


#### **APPLICATION:**

Flow cytometry; This antibody can be used as a negative isotypic control. The concentration will depend on condition.

Detailed procedure is provided in the following **PROTOCOLS**.

#### **INTENDED USE:**

For Research Use Only. Not for use in diagnostic procedures.



Flow cytometric analysis of mouse IgG2a reactivity on lymphocyte (R1), monocyte (R2) and granulocyte (R3).

# **PROTOCOLS:**

## Flow cytometric analysis for floating cells

We usually use Fisher tubes or equivalents as reaction tubes for all steps described below.

- 1) Wash the cells 3 times with washing buffer [PBS containing 2% fetal calf serum (FCS) and 0.1% NaN<sub>3</sub>].
- 2) Resuspend the cells with washing buffer  $(5x10^6 \text{ cells/mL})$ .
- 3) Add 50 μL of the cell suspension into each tube, and centrifuge at 500 x g for 1 minute at room temperature (20~25°C). Remove supernatant by careful aspiration.
- 4) Add 20 µL of Clear Back (human Fc receptor blocking reagent, MBL; code no. MTG-001) to the cell pellet after tapping. Mix well and incubate for 5 minutes at room temperature.
- 5) Add the isotype control antibody at the concentrations comparable to those of the specific antibody of interest. Mix well and incubate for 30 minutes at room temperature.
- 6) Add 1 mL of the washing buffer followed by centrifugation at 500 x g for 1 minute at room temperature. Remove supernatant by careful aspiration.
- 7) Resuspend the cells with 500  $\mu L$  of the washing buffer and analyze by a flow cytometer.

## Flow cytometric analysis for whole blood cells

We usually use Falcon tubes or equivalents as reaction tubes for all steps described below.

- 1) Add the isotype control antibody into each tube at the concentrations comparable to those of the specific antibody of interest.
- 2) Add 100  $\mu$ L of whole blood into each tube. Mix well and incubate for 30 minutes at room temperature (20~25 °C).
- 3) Add 1 mL of the washing buffer [PBS containing 2% fetal calf serum (FCS) and 0.1% NaN<sub>3</sub>] followed by centrifugation at 500 x g for 1 minute at room temperature. Remove supernatant by careful aspiration.
- 4) Lyse with OptiLyse C (for analysis on Beckman Coulter instruments) or OptiLyse B (for analysis on BD instruments), using the procedure recommended in the respective package inserts.
- 5) Add 1 mL of H<sub>2</sub>O to each tube and incubate for 10 minutes at room temperature.
- 6) Centrifuge at 500 x g for 1 minute at room temperature.
- 7) Add 1 mL of washing buffer followed by centrifugation at 500 x g for 1 minute at room temperature. Remove supernatant by careful aspiration.
- 8) Resuspend the cells with 500  $\mu L$  of the washing buffer and analyze by a flow cytometer.

# **RELATED PRODUCTS:**

[Functional grade antibody] M075-3M2 Mouse IgG1 isotype control FG (2E12) M076-3M2 Mouse IgG2a isotype control FG (6H3) M077-3M2 Mouse IgG2b isotype control FG (3D12) M080-3M2 Rat IgG1 isotype control FG (1H5) M081-3M2 Rat IgG2a isotype control FG (2H3) M090-3M2 Rat IgG2b isotype control FG (3G8) [Purified antibody] M075-3 Mouse IgG1 isotype control (2E12) Mouse IgG1 isotype control-FITC (2E12) M075-4 Mouse IgG1 isotype control-PE (2E12) M075-5 M075-8 Mouse IgG1 isotype control-Agarose (2E12) M075-A48 Mouse IgG1 isotype control-Alexa Fluor<sup>®</sup> 488 (2E12) M075-A64 Mouse IgG1 isotype control-Alexa Fluor® 647 (2E12) M076-3 Mouse IgG2a isotype control (6H3) M076-4 Mouse IgG2a isotype control-FITC (6H3) Mouse IgG2a isotype control-PE (6H3) M076-5 Mouse IgG2a isotype control-Alexa Fluor<sup>®</sup> 488 (6H3) M076-A48 Mouse IgG2a isotype control-Alexa Fluor<sup>®</sup> 647 (6H3) M076-A64 M077-3 Mouse IgG2b isotype control (3D12) M077-4 Mouse IgG2b isotype control-FITC (3D12) M077-5 Mouse IgG2b isotype control-PE (3D12) Mouse IgG2b isotype control-Alexa Fluor<sup>®</sup> 488 (3D12) M077-A48 Mouse IgG2b isotype control-Alexa Fluor® 647 (3D12) M077-A64 M078-3 Mouse IgG3 isotype control (6A3) M078-4 Mouse IgG3 isotype control-FITC (6A3) M079-3 Mouse IgM isotype control (7E10) M080-3 Rat IgG1 isotype control (1H5) M080-4 Rat IgG1 isotype control-FITC (1H5) M080-5 Rat IgG1 isotype control-PE (1H5) M080-A48 Rat IgG1 isotype control-Alexa Fluor<sup>®</sup> 488 (1H5) M081-3 Rat IgG2a isotype control (2H3) M081-4 Rat IgG2a isotype control-FITC (2H3) M081-5 Rat IgG2a isotype control-PE (2H3) M081-8 Rat IgG2a isotype control-Agarose (2H3) Rat IgG2a isotype control-Alexa Fluor® 488 (2H3) M081-A48 M090-3 Rat IgG2b isotype control (3G8) M090-4 Rat IgG2b isotype control-FITC (3G8) M090-5 Rat IgG2b isotype control-PE (3G8) Rat IgG2b isotype control-Alexa Fluor® 488 (3G8) M090-A48 M082-3 Rat IgG2c isotype control (6E12) M082-4 Rat IgG2c isotype control-FITC (6E12)

Normal Rabbit IgG-Agarose (polyclonal)

100510-2.1

#### LABEL LICENSES:

PM035-8

This product is provided under an agreement between LIFE TECHNOLOGIES Corporation, and Medical & Biological Laboratories Co., LTD, and the manufacture, use, sale or import of this product may be subject to one or more U.S. patents, pending applications, and corresponding non-U.S. equivalents, owned by or licensed to LIFE TECHNOLOGIES Corporation. For information on purchasing a license to this product for any other use, contact LIFE TECHNOLOGIES Corporation, Molecular Probes Detection Technologies, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, USA, Tel: (541) 465-8300. Fax: (541) 335-0504.