

Mouse anti-Human CD21

Clone NZ1

Cat. no. MON 1028-1

Specificity

MON1028 reacts with CD21, a 140 kD glycoprotein. The antigen is present on B-cells in sec. follicles of lymphnodes, tonsils and spleen; some malignant B-cells and on EBV trans-fected B-cell lines. This antibody does not block the biological function of CD21.

Immunoglobulin type

Mouse IgG2a_k

Use

The antibody can be used for immunostaining of cell suspensions (also after alcohol-acetic acid treatment), for immunohistochemistry of frozen sections and for immunoprecipitation.
Recommended dilution: 1:10-1:20 in PBS.

Presentation

1 ml purified ascites with approximately 1 mg IgG/ml. Add preservative if desired.

References

-

Storage

Store antibody at 2-8°C until expiry date.

FOR RESEARCH USE ONLY, NOT FOR DRUG, DIAGNOSTIC OR OTHER USE.