

Rat anti-Mouse IgG2a, clone LO-MG2a-3 (Monoclonal)

Clone no. LO-MG2a-3

MONOSAN

Product name	Rat anti-Mouse IgG2a, clone LO-MG2a-3 (Monoclonal)
Host	Rat
Applications	IHC-fr, ELISA
Species reactivity	mouse
Conjugate	-
Immunogen	mouse IgG1 heavy chain
Isotype	rat IgG2a k
Clonality	Monoclonal
Clone number	LO-MG2a-3
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	PBS with 0.1% sodium azide
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON24002

MONOSAN[®]

Rat anti-Mouse IgG2a, clone LO-MG2a-3 (Monoclonal)

Clone no. LO-MG2a-3

MONOSAN

Additional info

-

References

1. Denis O et al. Int Immunol 1993; 5(1):71-8
2. Pierre P et al. Euro J Immunol 1992; 22(12): 3179-82
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES