

Mouse anti-Human IgA, clone BFG-267 (Monoclonal)

Clone no. BFG-267

MONOSAN

Product name	Mouse anti-Human IgA, clone BFG-267 (Monoclonal)
Host	Mouse
Applications	IHC (1:40), ELISA (1:1000)
Species reactivity	human
Conjugate	-
Immunogen	IgA
Isotype	IgG1
Clonality	Monoclonal
Clone number	BFG-267
Size	1 ml
Concentration	0.5 mg
Format	-
Storage buffer	Purified and lyophilized
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human IgA, clone BFG-267 (Monoclonal)

Clone no. BFG-267

MONOSAN

Additional info

Purified and lyophilized from 0.01M sodium phosphate, 0.01M sodium borate, 0.15M sodium chloride, 1% dextran, 1%mannitol, pH 8. For Reconstitution use 1 ml sterile pure water. Specific against human IgA, a-chain. No reactivity to human IgE, IgG, IgM or other human plasma proteins.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES