

Mouse anti-Human IgA, clone BFG-267, HRP (Monoclonal)

Clone no. BFG-267

MONOSAN

---

Product name	Mouse anti-Human IgA, clone BFG-267, HRP (Monoclonal)
Host	Mouse
Applications	WB, ELISA, IHC
Species reactivity	human
Conjugate	HRP
Immunogen	IgA
Isotype	IgG1
Clonality	Monoclonal
Clone number	BFG-267
Size	1 ml
Concentration	1.2 mg/ml
Format	-
Storage buffer	0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human IgA, clone BFG-267, HRP (Monoclonal)

Clone no. BFG-267

MONOSAN

**Additional info**

Solution containing antibody 1.2 mg/ml in buffer (0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA). Specific against human IgA, a-chain. No reactivity to human IgE, IgG, IgM or other human plasma proteins.

**References**

1. -
2. -
3. -
4. -
5. -

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**