

Mouse anti-Human IgG, clone EFE-565 (Monoclonal)

Clone no. EFE-565

MONOSAN

Product name	Mouse anti-Human IgG, clone EFE-565 (Monoclonal)
Host	Mouse
Applications	WB, ELISA, IHC
Species reactivity	human
Conjugate	-
Immunogen	IgG
Isotype	IgG1
Clonality	Monoclonal
Clone number	EFE-565
Size	1 ml
Concentration	0.5 mg
Format	-
Storage buffer	Purified and lyophilized from 0.01M sodium phosphate, 0.01M sodium borate, 0.15M sodium chloride, 1% dextran, 1%mannitol, pH 8. Reconstitute
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human IgG, clone EFE-565 (Monoclonal)

Clone no. EFE-565

MONOSAN

Additional info

Purified and lyophilized from 0.01M sodium phosphate, 0.01M sodium borate, 0.15M sodium chloride, 1% dextran, 1%mannitol, pH 8. For Reconstitution use 1 ml sterile pure water. Human immunoglobulin IgG (gamma chain). No cross-reaction to human IgA, IgE, IgM or other human plasma proteins

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES