

Product datasheet MON5014P

MONOSAN[®]

Mouse anti-Human IgG, clone EFE-565, HRP (Monoclonal)

Clone no. EFE-565

MONOSAN

Product name	Mouse anti-Human IgG, clone EFE-565, HRP (Monoclonal)
Host	Mouse
Applications	WB, ELISA, IHC
Species reactivity	human
Conjugate	HRP
Immunogen	IgG
Isotype	IgG1
Clonality	Monoclonal
Clone number	EFE-565
Size	1 ml
Concentration	1.2 mg/ml
Format	-
Storage buffer	0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human IgG, clone EFE-565, HRP (Monoclonal)

Clone no. EFE-565

MONOSAN

Additional info

Solution containing antibody 1.2 mg/ml in buffer (0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA). Specific against human IgG, gamma-chain. No reactivity to human IgA, IgE, IgM or other human plasma proteins.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES