

Mouse anti-Human IgM, clone ICL-931 (Monoclonal)

Clone no. ICL-931

MONOSAN

Product name	Mouse anti-Human IgM, clone ICL-931 (Monoclonal)
Host	Mouse
Applications	WB, ELISA, IHC
Species reactivity	human
Conjugate	-
Immunogen	IgM
Isotype	IgG1
Clonality	Monoclonal
Clone number	ICL-931
Size	1 ml
Concentration	0.5 mg
Format	-
Storage buffer	Purified and lyophilized
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human IgM, clone ICL-931 (Monoclonal)

Clone no. ICL-931

MONOSAN

Additional info

Lyophilized from solution containing 0.25 mg of antibody, 1% mannitol, 0.15 M sodium chloride, 10 mM sodium phosphate, pH 7.5. For Reconstitution use 0.5 ml sterile pure water. Human immunoglobulin IgM (m chain). Does not react with IgA, IgE, IgG and other plasma proteins

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES