

Mouse anti-Human IgM, clone ICL-931, FITC (Monoclonal)

Clone no. ICL-931

MONOSAN

Product name	Mouse anti-Human IgM, clone ICL-931, FITC (Monoclonal)
Host	Mouse
Applications	IHC, ELISA, FC, EM
Species reactivity	human
Conjugate	FITC
Immunogen	IgM
Isotype	IgG1
Clonality	Monoclonal
Clone number	ICL-931
Size	1 ml
Concentration	1.2 mg/ml
Format	-
Storage buffer	0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human IgM, clone ICL-931, FITC (Monoclonal)

Clone no. ICL-931

MONOSAN

Additional info

Antibody solution 1 mg/ml in the buffer (0.05M sodium phosphate pH 7.5, 0.1 M NaCl) containing 50% glycerol and 5 mg/ml bovine serum albumin. Human immunoglobulin IgM. No cross-reactivity to human IgG, IgA and other human plasma proteins

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES