

Product datasheet PS100

MONOSAN[®]

Polyclonal Goat anti-Human IgG, affinity purified

Clone no. -

MONOSAN

Product name	Polyclonal Goat anti-Human IgG, affinity purified
Host	Goat
Applications	WB, ELISA (1:1000), IHC
Species reactivity	human
Conjugate	-
Immunogen	IgG
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol,pH 7.5
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Goat anti-Human IgG, affinity purified

Clone no. -

MONOSAN

Additional info

Solution containing antibody 1 mg/ml in buffer (0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, pH 7.5). Human immunoglobulin IgG. No reactivity against human IgA, IgM, other serum proteins.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES