

Product datasheet PS102F

MONOSAN[®]

Polyclonal Goat anti-Human IgA, FITC

Clone no. -

MONOSAN

Product name	Polyclonal Goat anti-Human IgA, FITC
Host	Goat
Applications	IHC (1:200), WB, ELISA, IF, FACS, EM
Species reactivity	human
Conjugate	FITC
Immunogen	IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1.2 mg/ml
Format	-
Storage buffer	0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Goat anti-Human IgA, FITC

Clone no. -

MONOSAN

Additional info

Solution 1.2 mg/ml in the buffer (0.05M sodium phosphate pH 7.5, 0.1 M NaCl) containing 50% glycerol and 5 mg/ml bovine. Human immunoglobulin IgA. No reactivity against human IgG, IgM, other serum proteins

References

1. Hecker J et al. Allergy 2012 Dec;67(12):1530-7
2. Hecker J et al. Mol Immunol. 2011 May;48(9-10):1236-44.
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES