

Product datasheet PS103F

MONOSAN[®]

Polyclonal Goat anti-Human IgG, IgM and IgA, FITC

Clone no. -

MONOSAN

| | |
|---------------------------|--|
| Product name | Polyclonal Goat anti-Human IgG, IgM and IgA, FITC |
| Host | Goat |
| Applications | WB, ELISA, IHC |
| Species reactivity | human |
| Conjugate | FITC |
| Immunogen | IgG, IgM and IgA |
| Isotype | - |
| Clonality | Polyclonal |
| Clone number | - |
| Size | 1 ml |
| Concentration | 1.2 mg/ml |
| Format | - |
| Storage buffer | 0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA |
| Storage until expiry date | -20°C |

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Goat anti-Human IgG, IgM and IgA, FITC

Clone no. -

MONOSAN

Additional info

Solution 1.2 mg/ml in the buffer (0.05M sodium phosphate pH 7.5, 0.1 M NaCl) containing 50% glycerol and 5 mg/ml bovine. Human immunoglobulins IgG, IgM, IgA. No reactivity against other serum proteins

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES