

Product datasheet PS104P

**MONOSAN**<sup>®</sup>

Polyclonal Rabbit anti-Human IgG, HRP

Clone no. -

MONOSAN

---

Product name	Polyclonal Rabbit anti-Human IgG, HRP
Host	Rabbit
Applications	WB, ELISA (1:1000), IHC
Species reactivity	human
Conjugate	HRP
Immunogen	IgG
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1.2 mg/ml
Format	-
Storage buffer	0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

PS104P

**MONOSAN**<sup>®</sup>

Polyclonal Rabbit anti-Human IgG, HRP

Clone no. -

MONOSAN

---

**Additional info**

Solution 1.2 mg/ml in the buffer (0.05M sodium phosphate pH 7.5, 0.1 M NaCl) containing 50% glycerol and 5 mg/ml bovine. Human immunoglobulin IgG. No cross-reactivity to human IgM, IgA and other human plasma proteins

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES