

Product datasheet PS106

**MONOSAN**<sup>®</sup>

Polyclonal Rabbit anti-Human IgA

Clone no. -

MONOSAN

---

Product name	Polyclonal Rabbit anti-Human IgA
Host	Rabbit
Applications	WB, ELISA (1:1000), IHC
Species reactivity	human, mouse (minor)
Conjugate	-
Immunogen	IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	50% glycerol, 0.1M NaCl, 0.01M sodium phosphate, pH 7.5
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Human IgA

Clone no. -

MONOSAN

**Additional info**

Solution 1mg/mL antibody, 50% glycerol, 0.1M NaCl, 0.01M sodium phosphate, pH 7.5. Human immunoglobulin IgA. No reactivity against human IgG, IgM, other serum proteins. It does have a minor cross-reactivity to mouse samples.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES