

Product datasheet PS106F

MONOSAN[®]

Polyclonal Rabbit anti-Human IgA, FITC

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-Human IgA, FITC
Host	Rabbit
Applications	WB, IHC (1:200), WB, IF, EM
Species reactivity	human, mouse (minor)
Conjugate	FITC
Immunogen	IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1.2 mg/ml
Format	-
Storage buffer	0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Human IgA, FITC

Clone no. -

MONOSAN

Additional info

Solution 1.2 mg/ml in the buffer (0.05M sodium phosphate pH 7.5, 0.1M NaCl, 50% glycerol, 5 mg/ml BSA). Human immunoglobulin IgA. No cross-reactivity to human IgG, IgM or other human plasma proteins. It does have a minor cross-reactivity to mouse samples.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES