

Product datasheet PS107

MONOSAN[®]

Polyclonal Rabbit anti-Human IgG, IgM and IgA

Clone no. -

MONOSAN

| | |
|---------------------------|---|
| Product name | Polyclonal Rabbit anti-Human IgG, IgM and IgA |
| Host | Rabbit |
| Applications | WB, ELISA, IHC |
| Species reactivity | human |
| Conjugate | - |
| Immunogen | IgG, IgM, IgA |
| Isotype | - |
| Clonality | Polyclonal |
| Clone number | - |
| Size | 1 ml |
| Concentration | 1 mg/ml |
| Format | - |
| Storage buffer | 50% glycerol, 0.1M NaCl, 0.01M sodium phosphate, pH 7.5 |
| Storage until expiry date | -20°C |

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

PS107

MONOSAN[®]

Polyclonal Rabbit anti-Human IgG, IgM and IgA

Clone no. -

MONOSAN

Additional info

Solution 1mg/mL antibody, 50% glycerol, 0.1M NaCl, 0.01M sodium phosphate, pH 7.5. Human immunoglobulins IgG, IgM, IgA. No reactivity against other serum proteins

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES