

Product datasheet PS108F

**MONOSAN**<sup>®</sup>

Polyclonal Rabbit anti-Goat IgG, IgM and IgA, FITC

Clone no. -

MONOSAN

---

Product name	Polyclonal Rabbit anti-Goat IgG, IgM and IgA, FITC
Host	Rabbit
Applications	WB, ELISA, IHC (1:40), IF, FACS, EM
Species reactivity	goat
Conjugate	FITC
Immunogen	IgG, IgM, IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1.2 mg/ml
Format	-
Storage buffer	0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Goat IgG, IgM and IgA, FITC

Clone no. -

MONOSAN

**Additional info**

Solution 1.2 mg/ml in the buffer (0.05M sodium phosphate pH 7.5, 0.1 M NaCl) containing 50% glycerol and 5 mg/ml bovine serum albumin. Binds goat immunoglobulins IgG, IgM, IgA. No cross reactivity with other goat serum proteins

**References**

1. -
2. -
3. -
4. -
5. -

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**