

Product datasheet PS108P

MONOSAN[®]

Polyclonal Rabbit anti-Goat IgG, IgM and IgA, HRP

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-Goat IgG, IgM and IgA, HRP
Host	Rabbit
Applications	WB, ELISA (1:10000), IHC
Species reactivity	goat
Conjugate	HRP
Immunogen	IgG, IgM, IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1.2 mg/ml
Format	-
Storage buffer	0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Goat IgG, IgM and IgA, HRP

Clone no. -

MONOSAN

Additional info

Liquid formulation containing antibody (1.2 mg/ml) in buffer (0.1 M NaCl, 0.05 M sodium phosphate, pH 7.5, 50% glycerol, 5 mg/ml BSA). Goat immunoglobulins IgG, IgM, IgA. Does not cross-react to other goat plasma proteins

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES