

Product datasheet PS110

**MONOSAN**<sup>®</sup>

Polyclonal Rabbit anti-Rat IgG, IgM and IgA

Clone no. -

MONOSAN

---

Product name	Polyclonal Rabbit anti-Rat IgG, IgM and IgA
Host	Rabbit
Applications	WB, ELISA (1:1000), IHC, DB (1:100)
Species reactivity	rat
Conjugate	-
Immunogen	IgG, IgM, IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	50% glycerol, 0.1 M NaCl, 50 mM sodium phosphate, pH 7.5
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Rat IgG, IgM and IgA

Clone no. -

MONOSAN

**Additional info**

Liquid formulation containing antibody (1 mg/ml) in 50% glycerol, 0.1 M NaCl, 50 mM sodium phosphate pH 7.5. Rat IgG, IgM and IgA. Does not cross-react with other rat plasma proteins

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES