

Product datasheet PS112F

**MONOSAN**<sup>®</sup>

Polyclonal Rabbit anti-Sheep IgG, IgM and IgA, FITC

Clone no. -

MONOSAN

---

Product name	Polyclonal Rabbit anti-Sheep IgG, IgM and IgA, FITC
Host	Rabbit
Applications	WB, ELISA, IHC (1:40), IF, FACS, EM
Species reactivity	sheep
Conjugate	FITC
Immunogen	IgG, IgM, IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1.2 mg/ml
Format	-
Storage buffer	0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Sheep IgG, IgM and IgA, FITC

Clone no. -

MONOSAN

**Additional info**

Solution 1.2 mg/ml in the buffer (0.05 M sodium phosphate pH 7.5, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA). Sheep immunoglobulins IgG, IgM, IgA. No cross-reactivity to other sheep plasma proteins.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES