

Product datasheet PS114

**MONOSAN**<sup>®</sup>

Polyclonal Rabbit anti-Horse IgG, IgM and IgA

Clone no. -

MONOSAN

---

Product name	Polyclonal Rabbit anti-Horse IgG, IgM and IgA
Host	Rabbit
Applications	WB, ELISA, IHC
Species reactivity	horse
Conjugate	-
Immunogen	IgG, IgM, IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	50% glycerol, 0.1 M NaCl, 0.05M sodium phosphate, pH 7.5
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

PS114

**MONOSAN**<sup>®</sup>

Polyclonal Rabbit anti-Horse IgG, IgM and IgA

Clone no. -

MONOSAN

---

**Additional info**

Solution 1 mg/mL antibody, 50% glycerol, 0.1M NaCl, 0.05M sodium phosphate, pH 7.5. Horse IgG, IgM, IgA

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES