

Product datasheet PS115

MONOSAN[®]

Polyclonal Rabbit anti-Chicken IgY

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-Chicken IgY
Host	Rabbit
Applications	WB, ELISA, IHC
Species reactivity	chicken
Conjugate	-
Immunogen	IgY
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	50% glycerol, 0.1 M NaCl, 50 mM sodium phosphate, pH 7.5
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Chicken IgY

Clone no. -

MONOSAN

Additional info

Liquid formulation containing antibody (1 mg/ml) in 50% glycerol, 50 mM sodium phosphate, 0.1M NaCl, pH 7.5. Chicken immunoglobulin IgY

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES