

Product datasheet PS115F

MONOSAN[®]

Polyclonal Rabbit anti-Chicken IgY, FITC

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-Chicken IgY, FITC
Host	Rabbit
Applications	WB, ELISA, IHC (1:200), IF, FACS, EM
Species reactivity	chicken
Conjugate	FITC
Immunogen	IgY
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1.2 mg/ml
Format	-
Storage buffer	0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Chicken IgY, FITC

Clone no. -

MONOSAN

Additional info

Solution 1.2 mg/ml in the buffer (0.1 M sodium chloride, 0.05M sodium phosphate pH 7.5, 50% glycerol, 5 mg/ml BSA). Binds chicken immunoglobulin IgY. Affinity to human IgM, IgG, IgA as measured by ELISA below 3%

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES