

Product datasheet PS117

**MONOSAN**<sup>®</sup>

Polyclonal Goat anti-Mouse IgG, IgM, IgA

Clone no. -

MONOSAN

---

Product name	Polyclonal Goat anti-Mouse IgG, IgM, IgA
Host	Goat
Applications	WB, ELISA (1:30000), IHC
Species reactivity	mouse
Conjugate	-
Immunogen	IgG, IgM, IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	50 mM sodium phosphate pH 7.4, 0.1 M NaCl, 50% glycerol
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Goat anti-Mouse IgG, IgM, IgA

Clone no. -

MONOSAN

**Additional info**

Liquid formulation containing antibody (1 mg/ml) in 50 mM sodium phosphate pH 7.4, 0.1 M NaCl, 50% glycerol. Affinity measured in ELISA: Mouse IgG 100%, IgM 90%, IgA 85%, murine plasma albumin <0.5%, human IgG, IgM, IgA <5%.

**References**

1. -
2. -
3. -
4. -
5. -

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**