

Polyclonal Goat anti-Mouse IgG, IgM, IgA - FITC

Clone no. -

MONOSAN

Product name	Polyclonal Goat anti-Mouse IgG, IgM, IgA - FITC
Host	Goat
Applications	WB, ELISA, IHC (1:200), IF, FACS, EM
Species reactivity	mouse
Conjugate	FITC
Immunogen	IgG, IgM, IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1.2 mg/ml
Format	-
Storage buffer	0.01 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Goat anti-Mouse IgG, IgM, IgA - FITC

Clone no. -

MONOSAN

Additional info

Solution 1.2 mg/ml in the buffer (0.01M sodium phosphate pH 7.5, 0.1M NaCl, 50% glycerol, 5 mg/ml BSA). Binds mouse (murine) immunoglobulins IgG, IgM, IgA. Affinity to human IgM, IgG, IgA as measured by ELISA below 3%

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES