

Product datasheet PS118

MONOSAN[®]

Polyclonal Goat anti-Rabbit IgG, IgM, IgA

Clone no. -

MONOSAN

Product name	Polyclonal Goat anti-Rabbit IgG, IgM, IgA
Host	Goat
Applications	WB, ELISA (1:1000), IHC, DB (1:100)
Species reactivity	rabbit
Conjugate	-
Immunogen	IgG, IgM, IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	50 mM sodium phosphate (pH 7.5), 100 mM NaCl, 50% glycerol
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

PS118

MONOSAN[®]

Polyclonal Goat anti-Rabbit IgG, IgM, IgA

Clone no. -

MONOSAN

Additional info

Supplied in 50 mM sodium phosphate (pH 7.5), 100 mM NaCl, 50% glycerol.
Rabbit IgG, IgM, IgA. Does not cross-react with other rabbit plasma proteins

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES