

Product datasheet PS121

MONOSAN[®]

Polyclonal Sheep anti-Rabbit IgG, IgM, IgA

Clone no. -

MONOSAN

Product name	Polyclonal Sheep anti-Rabbit IgG, IgM, IgA
Host	Sheep
Applications	WB, ELISA (1:1000), IHC, DB (1:100)
Species reactivity	rabbit
Conjugate	-
Immunogen	IgG, IgM, IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	10 mM sodium phosphate, 0.1 M NaCl, 50% glycerol
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

PS121

MONOSAN[®]

Polyclonal Sheep anti-Rabbit IgG, IgM, IgA

Clone no. -

MONOSAN

Additional info

Liquid formulation containing antibody (1 mg/ml) in 50%, 0.1 M NaCl, 10 mM sodium phosphate, pH 7.5

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES