

Product datasheet PS122

**MONOSAN**<sup>®</sup>

Polyclonal Sheep anti-Human IgG, IgM, IgA

Clone no. -

MONOSAN

---

Product name	Polyclonal Sheep anti-Human IgG, IgM, IgA
Host	Sheep
Applications	WB, ELISA (1:1000), IHC
Species reactivity	human
Conjugate	-
Immunogen	IgG, IgM, IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	0.01 M sodium phosphate, 0.1 M NaCl, 50% glycerol
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

PS122

**MONOSAN**<sup>®</sup>

Polyclonal Sheep anti-Human IgG, IgM, IgA

Clone no. -

MONOSAN

---

**Additional info**

Solution 1 mg/ml in the buffer (0.01M sodium phosphate pH 7.5, 0.1M NaCl, 50% glycerol). Human immunoglobulins IgG, IgM, IgA

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES