

Product datasheet PS124

MONOSAN[®]

Polyclonal Sheep anti-Human IgA

Clone no. -

MONOSAN

Product name	Polyclonal Sheep anti-Human IgA
Host	Sheep
Applications	WB, ELISA (1:10000), IHC
Species reactivity	human
Conjugate	-
Immunogen	IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	0.05 M sodium phosphate, 50% glycerol, pH 7.5
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

PS124

MONOSAN[®]

Polyclonal Sheep anti-Human IgA

Clone no. -

MONOSAN

Additional info

Aqueous solution containing the antibody 1 mg/ml, 50% glycerol, 0.05 M sodium phosphate pH 7.5. Human IgA – 100%, human IgG <10%, human IgM <5%, human serum albumin <1%

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES