

Product datasheet PS127

MONOSAN[®]

Polyclonal Sheep anti-Human IgM

Clone no. -

MONOSAN

Product name	Polyclonal Sheep anti-Human IgM
Host	Sheep
Applications	IHC, ELISA (1:30000), WB
Species reactivity	human
Conjugate	-
Immunogen	IgM
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	0.5 ml
Concentration	2 mg/ml
Format	-
Storage buffer	50% glycerol, 50 mM sodium phosphate, pH 7.5
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Sheep anti-Human IgM

Clone no. -

MONOSAN

Additional info

Liquid formulation containing antibody (2 mg/ml) in 50% glycerol, 50 mM sodium phosphate, pH 7.5. Affinity measured in ELISA: human IgM 100%, human IgA <2%, human IgG <1%, human albumin <0.2%, murine IgG, IgM, IgA <1%.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES