

Mouse anti-Human Melanoma Associated Antigen

Clone PNL2

Cat. no. MON 3307

Specificity

MON 3307 is useful as an immunohistochemical reagent to stain melanocytes and tumors derived therefrom. Anti-PNL2 reactivity is identified in the cytoplasm of cutaneous and oral mucosal melanocytes. Anti-PNL2 labels intraepidermal nevi, while the dermal components of compound nevi are largely non-reactive.

Immunoglobulin type

Mouse IgG1

Use

This antibody can be used in IHC Paraffin sections (1:25 – 1:100)

Presentation

1 ml supernatant in Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide

References

- Rochaix P, et al. Mod Pathol. 2003; 16:481-90.
- Klaus J, et al. Am J Surg Pathol. 2005; 29:400-6.
- Morris LG, et al. Head Neck. 2008; 30:771-5.
- Zhe X, et al. J Histo and Cyto. 2004; 52:1537-42

Storage

Store antibody at 2-8°C until expiry date.

FOR RESEARCH USE ONLY, NOT FOR DRUG, DIAGNOSTIC OR OTHER USE.

11/19 MONOSAN Cat.# MON 3307

FOR-044 23-10-2019